# Testera Multi-Drug 3 Drugs Rapid Test (Oral Fluid)

## Package Insert

REF DSD-833-D English

### Instruction Sheet for testing of any combination of the following drugs:

THC/COT/K2

A rapid test for the simultaneous, qualitative detection of multiple drugs or drug metabolites in human oral fluid. For healthcare professionals including professionals at point of care sites. Immunoassay for in vitro diagnostic use only.

#### 【INTENDED USE】

The Multi-Drug Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of multiple drugs or drug metabolites in human oral fluid at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	50
Cotinine (COT)	Cotinine	20
Synthetic Marijuana (K2)	JWH-018, JWH-073	25

This assay provides only a preliminary test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) is the preferred confirmatory methods. Professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

#### 【SUMMARY】

The Multi-Drug Rapid Test is a rapid, oral fluid screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in human oral fluid.

#### Marijuana (THC)

11-nor- $\Delta^9$ -tetrahydrocannabinol-9-carboxylic acid ( $\Delta^9$ -THC-COOH), the metabolite of THC ( $\Delta^9$ -tetrahydrocannabinol), is detectable in oral fluid shortly after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of the drug in the buccal cavity.<sup>1</sup> Historical studies have shown a window of detection for THC in oral fluid of up to 14 hours after drug use.<sup>1</sup>

#### Cotinine (COT)

Cotinine is the first-stage metabolite of nicotine, a toxic alkaloid that produces stimulation of the autonomic ganglia and central nervous system when in humans. Nicotine is a drug to which virtually every member of a tobacco-smoking society is exposed whether through direct contact or second-hand inhalation. In addition to tobacco, nicotine is also commercially available as the active ingredient in smoking replacement therapies such as nicotine gum, transdermal patches and nasal sprays.

Although nicotine is excreted in oral fluid, the relatively short half-life of the drug makes it an unreliable maker for tobacco use. Cotinine, however, demonstrates a substantially longer half-life than nicotine bears a high correlation with plasma cotinine levels and has been found to be the best maker for smoking status compared with oral fluid nicotine measurement, breath carbon monoxide testing and plasma thiocyanate testing. The window of detection for cotinine in oral fluid test is expected to be up to 1-2 days after nicotine use.

#### Synthetic Marijuana (K2)

Synthetic Marijuana or K2 is a psychoactive herbal and chemical product that, when consumed, mimics the effects of Marijuana. It is best known by the brand names K2 and Spice, both of which have largely become genericized trademarks used to refer to any synthetic Marijuana product. The studies suggest that synthetic marijuana intoxication is associated with acute psychosis, worsening of previously stable psychotic disorders, and also may have the ability to trigger a chronic (long-term) psychotic disorder among vulnerable individuals such as those with a family history of mental illness.<sup>2</sup>

Elevated levels of oral fluid metabolites are found within hours of exposure and remain detectable window up to 24-48 hours after smoking (depending on usage/dosage).

#### 【ASSAY PRINCIPLE】

The Multi-Drug Rapid Test is an immunoassay based on the principle of competitive binding. Drugs that may be present in the oral fluid specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a portion of the oral fluid specimen migrates upward by capillary action. A drug, if present in the oral fluid specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region of the specific drug strip. The presence of drug above the cut-off concentration in the oral fluid specimen will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region.

A drug-positive oral fluid specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative oral fluid specimen will generate

a line in the test line region because of the absence of drug competition.

To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

#### 【REAGENTS】

Each test contains membrane strips coated with drug-protein conjugates (purified bovine albumin) on the test line, a goat polyclonal antibody against gold-protein conjugate at the control line, and a dye pad which contains colloidal gold particles coated with mouse monoclonal antibody specific to corresponding drug.

#### 【PRECAUTIONS】

- Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- Oral fluid is not classified as biological hazard unless derived from a dental procedure.
- The used Device should be discarded according to local regulations.

#### 【STORAGE AND STABILITY】

Store as packaged in the sealed pouch at 2-30 °C. The test is stable through the expiration date printed on the sealed pouch. The test Device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

#### 【SPECIMEN COLLECTION AND PREPARATION】

The oral fluid specimen should be collected using the device with the kit. Follow the detailed Directions for Use below. No other collection device should be used with this assay. Oral fluid collected at any time of the day may be used.

#### 【MATERIALS】

- Materials Provided**
- Test Devices
  - Package insert
- Materials Required but Not Provided**

- Timer

#### 【DIRECTIONS FOR USE】

**Allow the test device, specimen and/or controls to reach room temperature (15-30 °C) prior to testing. Instruct the donor to not place anything in the mouth including food, drink, gum or tobacco products for at least 10 minutes prior to collection.**

- Bring the pouch to room temperature before opening it. Remove the test from the sealed pouch and use it within one hour.
- Take off the Device cap and collect oral fluid specimen as follows.

**Important:** Place the absorbent wick against the upper, lower jaws and roots to enrich the oral fluid. Insert the sponge end into the mouth, actively swab around the gums on both sides of the mouth (10-15 times) to assist saturation.


Put the absorbent wick under the tongue to collect oral fluid until the flow appear in the test windows (approximately 60 seconds) and then take out the device and start a timer.

If no flow appeared repeat the procedure in steps above until the flow appear. If no flow appeared after triplicate of steps above, discard the device, review procedures with the donor and repeat the test using a new device.

- Place the test device on a clean and level surface.

- Read the test result at **3-10 minutes**.

If all lines are clearly visible at 3 minutes or sooner, then the test can be interpreted as negative and discarded. **If any lines are not visible at 3 minutes**, then the test should be re-read at 10 minutes.

#### 【INTERPRETATION OF RESULTS】

(Please refer to the previous illustration)

**NEGATIVE:** A colored line appears in the control region (C) and another colored line appears in the test region (T). This negative result means that the concentrations in the oral fluid sample are below the designated cut-off levels for a particular drug tested.

**\*NOTE:** The shade of the colored lines(s) in the Test regions (T) may vary. The result should be considered negative whenever there is even a faint line.

**POSITIVE:** A colored line appears in the control region (C) and no line appears in the test region (T). The positive result means that the drug concentration in the oral fluid sample is greater than the designated cut-off for a specific drug.

**INVALID:** No line appears in the control region (C). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for Control line failure. Read the directions

again and repeat the test with a new test device. If the result is still invalid, contact your manufacturer.

#### 【QUALITY CONTROL】

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

#### 【LIMITATIONS】

- The Multi-Drug Rapid Test provides only a qualitative, preliminary result. A secondary analytical method must be used to obtain a confirmed result. Gas Chromatography/Mass Spectrometry (GC/MS) is preferred confirmatory methods.<sup>3</sup>
- A positive test result does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate a drug-free specimen. Drug may be present in the specimen below the cutoff level of the assay.

#### 【PERFORMANCE CHARACTERISTICS】

##### Analytical Sensitivity

A Phosphate-buffered saline (PBS) pool was spiked with drugs to target concentrations of  $\pm 50\%$  cut-off,  $\pm 25\%$  cut-off and  $+300\%$  cut-off and tested with the Multi-Drug Rapid Test. The results are summarized below.

Drug Concentration Cut-off Range	THC 50		COT 20		K2 25	
	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0
-25% Cut-off	27	3	25	5	26	4
Cut-off	12	18	20	10	15	15
+25% Cut-off	8	22	7	23	3	27
+50% Cut-off	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30

##### Analytical Specificity

The following table lists the concentration of compounds (ng/mL) above which the Multi-Drug Rapid Test identified positive results at a read time of 3-10 minutes.

Compound	ng/mL	Compound	ng/mL
<b>Marijuana (THC 50)</b>			
11-nor- $\Delta^9$ -THC-9 COOH	50	$\Delta^8$ -THC	25,000
Cannabinol	50,000	$\Delta^9$ -THC	40,000
11-nor- $\Delta^8$ -THC-9 COOH	40		
<b>Cotinine (COT 20)</b>			
(-)-Cotinine	20	(-)-Nicotine	300
<b>Synthetic Marijuana (K2-25)</b>			
JWH-018 5-Pentanoic acid metabolite	25	JWH-018 4-Hydroxypentyl metabolite	200
JWH-073 4-butanoic acid metabolite	25	JWH-018 5-Hydroxypentyl metabolite	250
JWH-073 4-Hydroxybutyl metabolite	250		

##### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following compounds demonstrated no false positive results on the Multi-Drug Rapid Test when tested with at concentrations up to 100  $\mu\text{g/mL}$ .


Acetaminophen	d/l-Chlorpheniramine	Sulfamethazine
N-Acetylprocainamide	Chloroquine	Tetracycline
Aminopyrine	Clonidine	Tetrahydrocortisone $\beta$ (3-D-glucuronide)
Ampicillin	l-Cotinine	Thioridazine
Apomorphine	Deoxycorticosterone	Tolbutamide
Atropine	Diclofenac	Trifluoperazine
Benzoic acid	Digoxin	d/l-Tryptophan
d/l-Brompheniramine	l- $\Psi$ -Ephedrine	Uric acid
Chloral-hydrate	Estrone-3-sulfate	Ketoprofen
Chlorothiazide	l(-)-Epinephrine	Loperamide
Chlorpromazine	Fenoprofen	Meprobamate
Cholesterol	Gentisic acid	Nalidixic acid
Cortisone	Hydralazine	Niacinamide
Creatinine	Hydrocortisone	Norethindrone
Dextromethorphan	p-Hydroxytyramine	Noscapine
Diflunisal	Iproniazid	Oxalic acid
Diphenhydramine	Isosuprine	Oxymetazoline
$\beta$ -Estradiol	Labeltalol	Penicillin-G
Ethyl-p-aminobenzoate	Meperidine	Phenazazine
Erythromycin	Methylphenidate	Trans-2-phenylcyclopropylamine

Furosemide	Naproxen	hydrochloride
Hemoglobin	Nifedipine	Prednisolone
Hydrochlorothiazide	d-Norpropoxyphene	d/l-Propranolol
o-Hydroxyhippuric acid	d/l-Octopamine	d-Pseudoephedrine
Ibuprofen	Oxolinic acid	Quinine
d/l-Isoproterenol	Papaverine	Ranitidine
Acetophenetidin	Pentazocine hydrochloride	Serotonin
Acetylsalicylic acid	Phenelzine	Sulindac
Amoxicillin	Phenylpropanolamine	Tetrahydrocortisone 3-acetate
l-Ascorbic acid	Prednisone	Thiamine
Aspartame	d-Propoxyphene	d/l-Tyrosine
Benzilic acid	Quinacrine	Triamterene
Benzphetamine	Quindine	Trimethoprim
Caffeine	Salicylic acid	Tyramine
Chloramphenicol	Zomepirac	Verapamil

**【BIBLIOGRAPHY】**

- Schramm, W. et al, "Drugs of Abuse in Saliva: A Review," *J Anal Tox*, 1992 Jan-Feb; 16 (1), pp 1-9.
- Lane JC, Tennon MB, Lawless ST, et al. Movement disorder after withdrawal of fentanyl infusion. *J Pediatr*, 1991, 119(4): 649-651.
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

**Index of Symbols**

	Consult instructions for use		Contains sufficient for <n> test		Authorized representative in the European Community/European Union
	In vitro diagnostic medical device		Use-by date		Do not reuse
	Store between 2-30 °C		Batch code		Catalogue number
	Do not use if package is damaged and consult instructions for use		Manufacturer		

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MedNet EC-REP GmbH