Testera Multi-Drug 6 Drugs Rapid Test Panel (Urine)

Package Insert

REF DOA-164 English

Instruction Sheet for testing of any combination of the following drugs: COC/THC/MET/MDMA/FYL/MEP

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine. For healthcare professionals including professionals at point of care sites. Immunoassay for *in vitro* diagnostic use only.

INTENDED USE

The Multi-Drug Rapid Test Panel is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in human urine at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
Cocaine (COC)	Benzoyllecgonine	300
Marijuana (THC)	11-nor- Δ^8 -THC-9 COOH	50
Methamphetamine (MET)	d-Methamphetamine	1,000
Methylenedioxyamphetamine (MDMA)	d,l-Methylenedioxy-methamphetamine	500
Fentanyl (FYL)	Norfentanyl	20
Mephedrone (MEP)	Mephedrone	100

Configurations of the Multi-Drug Rapid Test Panel come with any combination of the above listed drug analytes. This assay provides only a preliminary test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

PRINCIPLE

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test region of the specific drug dipstick. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test region.

A drug-positive urine specimen will not generate a colored line in the specific test region of the dipstick because of drug competition, while a drug-negative urine specimen will generate a line in the test region because of the absence of drug competition.

To serve as a procedural control, a colored line will always appear at the control region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test line contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

PRECAUTIONS

- For healthcare professionals including professionals at point of care sites.
- Immunoassay for *in vitro* diagnostic use only. The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The Test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay
The urine specimen should be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing.


MATERIALS

- Test Panels
- Materials Provided
 - Package insert
- Materials Required But Not Provided
 - Specimen collection containers
- Timer

DIRECTIONS FOR USE

Allow the test, urine specimen and/or controls to reach room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test panel from the sealed pouch and use it within one hour.
- Remove the cap.
- With the arrow pointing toward the urine specimen, immerse the test panel vertically in the urine specimen for at least 10 to 15 seconds. **Immerse the dipstick to at least the level of the wavy lines, but not above the arrow on the test panel.**
- Replace the cap and place the test panel on a non-absorbent flat surface.
- Start the timer and wait for the colored line(s) to appear.
- The drug strip result should be read at **5 minutes**. Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration)

NEGATIVE: A colored line appears in the control region (C) and another colored line appears in the test region (T). This negative result means that the concentrations in the urine sample are below the designated cut-off levels for a particular drug tested.

NOTE: The shade of the colored lines(s) in the test region (T) may vary. The result should be considered negative whenever there is even a faint line.

POSITIVE: A colored line appears in the control region (C) and no line appears in the test region (T). The positive result means that the drug concentration in the urine sample is greater than the designated cut-off for a specific drug.

INVALID: No line appears in the control region (C). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Read the directions again and repeat the test with a new test. If the result is still invalid, contact your manufacturer.

QUALITY CONTROL

A procedural control is included in the test. A line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The Multi-Drug Rapid Test Panel provides only a qualitative, preliminary result. A secondary analytical method must be used to obtain a confirmed result. Gas Chromatography /Mass Spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may cause erroneous results.
- A positive result does not indicate level or intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- This test does not distinguish between drugs of abuse and certain medications.
- A positive test result may be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

	Accuracy					
	% Agreement with GC/MS					
	COC 300	THC 50	MET 1,000	MDMA 500	FYL 20	MEP 100
Positive Agreement	98.2%	97.9%	96.2%	98.1%	98.8%	90.5%
Negative Agreement	97.8%	98.1%	97.1%	99.3%	99.4%	97.0%
Total Results	98.0%	98.0%	96.8%	98.8%	99.2%	95.4%

	% Agreement with Commercial Kit					
	COC 300	THC 50	MET 500	MDMA 500	MEP 100	FYL 20
Positive Agreement	>99.9%	>99.9%	>99.9%	>99.9%	*	*
Negative Agreement	>99.9%	>99.9%	>99.9%	>99.9%	*	*
Total Results	>99.9%	>99.9%	>99.9%	>99.9%	*	*

*Note: Based on GC/MS data instead of Commercial Kit.

Precision

A study was conducted at three hospitals using three different lots of product to demonstrate the within run, between run and between operator precision. An identical card of coded specimens, containing drugs at concentrations of negative, 50% and 25% cut-off level, was labeled, blinded and tested at each site. **The results gained \cong 75% accuracy in \pm 25% cut-off level specimen and 100% accuracy in negative and \pm 50% cut-off level specimen.**

Analytical Sensitivity

A drug-free urine pool was spiked with drugs at the listed concentrations. The results are summarized below.

Drug Concentration Cut-off Range	COC 300		THC 50		MET 1,000		MDMA 500		FYL 20		MEP 100	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0

-50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	26	4	26	4	27	3	25	5	27	3	27	3
Cut-off	13	17	14	16	16	14	14	16	14	16	17	13
+25% Cut-off	3	27	3	27	3	27	4	26	4	26	5	25
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30
300% Cut-off	0	30	0	30	0	30	0	30	0	30	0	30

Analytical Specificity

The following table lists the concentrations of compounds (ng/mL) that are detected as positive in urine by the Multi-Drug Rapid Test at 5 minutes.

Analytes	conc. (ng/mL)	Analytes	conc. (ng/mL)
COCAINE (COC 300)			
Benzoyllecgonine	300	Cocaethylene	20,000
Cocaine HCl	200	Ecgonine	30,000
MARIJUANA (THC 50)			
Cannabinol	35,000	Δ^8 -THC	17,000
11-nor- Δ^8 -THC-9 COOH	30	Δ^9 -THC	17,000
11-nor- Δ^9 -THC-9 COOH	50		
METHAMPHETAMINE (MET 1,000)			
p-Hydroxymethamphetamine	25,000	(\pm)-3,4-Methylenedioxy-methamphetamine	12,500
D-Methamphetamine	1,000		
L-Methamphetamine	20,000	Mephentermine	50,000
METHYLENEDIOXYMETHAMPHETAMINE (MDMA 500) Ecstasy			
(\pm) 3,4-Methylenedioxy-methamphetamine HCl	500	3,4-Methylenedioxyethylamphetamine	300
(\pm) 3,4-Methylenedioxyamphetamine HCl	3,000		
FENTANYL (FYL 20)			
Alfentanyl	600,000	Bupirone	15,000
Fenfluramine	50,000	Fentanyl	100
Norfentanyl	20	Sufentanyl	50,000
paliperidone	1,250	Risperidone	5,000
MEPHEDRONE (MEP 100)			
Mephedrone HCl	100	R(+)-Methcathinone HCl	1,500
S(-)-Methcathinone HCl	500	3-Fluoromethcathinone HCl	1,500
4-Fluoromethcathinone HCl	300	Methoxyphenamine	100,000

Effect of Urinary Specific Gravity

Fifteen (15) urine samples of normal, high and low specific gravity ranges (1.005-1.045) were spiked with drugs at 50% below and 50% above cut-off levels respectively. The Multi-Drug Rapid Test was tested in duplicate using fifteen drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with drugs at 50% below and 50% above cut-off levels. The spiked, pH-adjusted urine was tested with the Multi-Drug Rapid Test. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug positive urine containing above related calibrator substances. The following compounds show no cross-reactivity when tested with the Multi-Drug Rapid Test at a concentration of 100 μ g/mL.


Non Cross-Reacting Compounds

Acetophenetidin	Cortisone	Zomepirac	Quinidine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinine
Acetylsalicylic acid	Deoxycorticosterone	Labetalol	Salicylic acid
Aminopyrine	Dextromethorphan	Loperamide	Serotonin
Amoxicillin	Diclofenac	Meprobamate	Sulfamethazine
Ampicillin	Diffunisal	Isoxsuprine	Sulindac
l-Ascorbic acid	Digoxin	d,l-Propranolol	Tetracycline
Apomorphine	Diphenhydramine	Nalidixic acid	Tetrahydrocortisone, 3-acetate
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone
Atropine	β -Estradiol	Niacinamide	Tetrahydrozoline
Benzilic acid	Estrone-3-sulfate	Nifedipine	Thiamine
Benzoic acid	Erythromycin	Norethindrone	Thioridazine
Bilirubin	Fenopropfen	Noscapine	d,l-Tyrosine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Tolbutamide
Cannabidiol	Genistic acid	Oxalic acid	Triamterene
Chloral hydrate	Hemoglobin	Oxolinic acid	Trifluoperazine
Chloramphenicol	Hydralazine	Oxymetazoline	Trimethoprim
Chlorothiazide	Hydrochlorothiazide	Penicillin-G	d,l-Tryptophan
d,l-Chlorpheniramine	Hydrocortisone	Perphenazine	Uric acid
Chlorpromazine	o-Hydroxyhippuric acid	Phenelzine	Verapamil
Cholesterol	3-Hydroxytyramine	Prednisone	
Clonidine	d,l-Isoproterenol		

BIBLIOGRAPHY

- Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 6th Ed. Biomedical Publ., Foster City, CA 2002.

Index of Symbols

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Temperature limit
IVD	<i>In vitro</i> diagnostic medical device	LOT	Batch code	REF	Catalogue number
EC REP	Authorized representative in the European Community/European Union		Use-by date		Do not re-use
	Do not use if package is damaged and consult instructions for use		Manufacturer		Caution

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MedNet EC-REP GmbH