Testera Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma)
Package Insert

REF IHVGM-402 English

A rapid test for the qualitative detection of IgG and IgM antibodies to Hantavirus in human whole blood, serum or plasma specimens. For professional in vitro diagnostic use only.

【INTENDED USE】

The Hantavirus IgG/IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of Hantavirus IgG/IgM in human whole blood, serum or plasma specimens.

【SUMMARY】

Hantaviruses are enveloped, single-stranded RNA viruses belonging to the family Bunyaviridae(genus Hantavirus) that cause two major clinical syndromes in humans:[1,3,5,6,7] hemorrhagic fever with renal syndrome (HFRS) in Europe and Asia, and hantavirus pulmonary syndrome (HPS), also known as hantavirus cardiopulmonary syndrome (HCPS), in the Americas.[1,2,4,7] These zoonotic pathogens are primarily transmitted to humans through inhalation of aerosolized excreta from persistently infected rodent reservoirs.[1,5] With case fatality rates ranging from 5–15% for HFRS to 20–40% for HPS, hantaviruses represent a significant public health concern in endemic regions.[1,2,3,4,7] The Andes virus is unique among hantaviruses in its documented capability for limited person-to-person transmission.[1,8]

【PRINCIPLE】

The Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is a qualitative, membrane based immunoassay for the detection of IgG and IgM antibodies to Hantavirus in whole blood, serum or plasma. The membrane is pre-coated with mouse anti-human IgG and mouse anti-human IgM on the test line region of the cassette. During testing, the whole blood, serum or plasma specimen reacts with recombinant Hantavirus antigen conjugated colloid gold. The mixture migrates upward on the membrane chromatographically by capillary action to react with mouse anti-human IgG and mouse anti-human IgM on the membrane and generate a colored line. Presence of this colored line indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【REAGENTS】

The test contains recombinant Hantavirus antigen conjugated colloid gold, mouse anti-human IgG and mouse anti-human IgM coated on the membrane.

【PRECAUTIONS】

1. For professional *in vitro* diagnostic use only. Do not use after expiration date.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Do not use test if pouch is damaged.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
5. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
6. The used test should be discarded according to local regulations.
7. Humidity and temperature can adversely affect results.

【STORAGE AND STABILITY】

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【SPECIMEN COLLECTION AND PREPARATION】

- The Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or finger prick), serum or plasma.
- To collect **Finger Prick Blood Specimens**:
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
 - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
 - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
 - Add the Finger Prick Blood specimen to the test cassette by using a dropper or micropipette measuring 20μL. The dropper provided with the test dispenses approximately 20μL in two drops even if more blood is aspirated in the dropper.
- If you plan to use serum or plasma specimens, please separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear non-hemolyzed specimens.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days, for long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by finger prick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiological agents.
- EDTA K2, Heparin sodium, Citrate sodium and Potassium Oxalate can be used as the anticoagulant for collecting the specimen.

【MATERIALS】

Materials provided: • Test Cassettes • Droppers • Buffer • Package insert

Materials required but not provided: Specimen Collection Containers • Centrifuge • Timer • Lancets (for Finger Prick Blood only) • Pipette

【DIRECTIONS FOR USE】

Allow the test, specimen, buffer and/or controls to reach room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface

For **Serum or Plasma** specimen:

• Use a dropper: Hold the dropper vertically and transfer **1 drop of specimen (approximately 10 μL)** to the specimen well (S) of the test cassette, then add **2 drops of buffer (approximately 80 μL)** into the buffer well (B) and start the timer.

Avoid trapping air bubbles in the specimen well (S).

• Use a micropipette: Pipette and dispense **10 μL of specimen** to the specimen well **(S)** of the test cassette, then add **2 drops of buffer (approximately 80 μL**) into the buffer well **(B)** and start the Timer.

For **Venipuncture and Finger Prick Blood** specimen:


• Use a dropper: Hold the dropper vertically and transfer **2 drops of specimen (approximately 20 μL)** to the specimen well **(S)** of the test cassette, then add **2 drops of buffer (approximately 80 μL)** into the buffer well **(B)** and start the timer.

Avoid trapping air bubbles in the specimen well (S).

• Use a micropipette: Pipette and dispense **20 μL of specimen** to the specimen well **(S)** of the test cassette, then add **2 drops of buffer (approximately 80 μL**) into the buffer well **(B)** and start the Timer.

3. Wait for the colored line(s) to appear. **Read the results at 10 minutes**, do not interpret the results after 20 minutes.

Note: It is suggested not to use the buffer beyond 6 months after opening the vial.

【INTERPRETATION OF RESULTS】

IgG POSITIVE: Two colored lines appear. One colored line should be in the control region (C) and another colored line should be in the IgG region.

IgM POSITIVE: Two colored lines appear. One colored line should be in the control region (C) and another colored line should be in the IgM region.

IgG and IgM POSITIVE: Three colored lines appear. One colored line should be in the control region (C) and another two colored lines should be in the IgG and IgM regions.

NOTE: The intensity of the color in the IgG and IgM region(s) will vary depending on the concentration of Hantavirus antibodies present in the specimen. Therefore, any shade of red in the IgG and IgM region(s) should be considered positive.

NEGATIVE: One colored line appears in the control region (C). No colored line appears in the IgG and IgM region(s).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

1. The Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is for *in vitro* diagnostic use only. This test should be used for qualitative detection of Hantavirus IgG/IgM antibodies in whole blood, serum or plasma specimens. Neither the quantitative value nor the rate of increase in the concentration of Hantavirus IgG/IgM antibodies can be determined by this qualitative test.
2. The Hantavirus IgG/IgM Rapid Test Cassette (Whole blood/Serum/Plasma) will only indicate the presence of Hantavirus IgG/IgM antibodies in the specimen and should not be used as the sole criterion for the diagnosis of Hantavirus infections.
3. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
4. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is suggested. A negative result at any time does not preclude the possibility of Hantavirus infection.
5. The hematocrit level of the whole blood can affect the test results. Hematocrit level needs to be between 25% and 65% for accurate results.

【PERFORMANCE CHARACTERISTICS】

Sensitivity and Specificity

The Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) was compared with a leading commercial ELISA Hantavirus IgG/IgM tests; the results show that Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) has a high sensitivity and specificity.

IgG Results

Rapid Test / ELISA	Positive	Negative	Total Results
Positive	58	3	61
Negative	2	162	164
Total Results	60	165	225

Relative Sensitivity: 96.7% (95%CI*: 88.5%-99.6%) ; Relative Specificity: 98.2% (95%CI*: 94.8%-99.6%); Accuracy: 97.8% (95%CI*:94.9%-99.3%) *Confidence Interval

IgM Results

Rapid Test / ELISA	Positive	Negative	Total Results
Positive	51	2	53
Negative	3	169	172
Total Results	54	171	225

Relative Sensitivity: 94.4% (95%CI*: 84.6%-98.8%) ;Relative Specificity: 98.8% (95%CI*: 95.8%-99.9%); Accuracy: 97.8% (95%CI*:94.9%-99.3%) *Confidence Interval

Precision

Intra-Assay. Within-run precision has been determined by using 10 replicates of five specimens: a negative, a Hantavirus IgM low titer positive, a Hantavirus IgM high titer positive, a Hantavirus IgG low titer positive and a Hantavirus IgG high titer positive. The negative, a Hantavirus IgM low titer positive, a Hantavirus IgM high titer positive, a Hantavirus IgG low titer positive and a Hantavirus IgG high titer positive values were correctly identified 100% of the time.

Inter-Assay. Between-run precision has been determined by 10 independent assays on the same five specimens: a negative, a Hantavirus IgM low titer positive, a Hantavirus IgM high titer positive, a Hantavirus IgG low titer positive and a Hantavirus IgG high titer positive. Three different lots of the Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) have been tested over a 3-days period using negative, a Hantavirus IgM low titer positive, a Hantavirus IgM high titer positive, a Hantavirus IgG low titer positive and a Hantavirus IgG high titer positive specimens. The specimens were correctly identified 100% of the time.

Cross-reactivity

The Hantavirus IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested with HAMA, RF, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, Syphilis, HIV, *H. Pylori*, MONO, CMV, Rubella and TOXO positive specimens. The results showed no cross-reactivity.

Interfering Substances

The following potentially interfering substances were added to Hantavirus negative and positive specimens.

Acetaminophen: 20 mg/dL	Caffeine: 20 mg/dL
Acetylsalicylic Acid: 20 mg/dL	Gentisic Acid: 20 mg/dL
Ascorbic Acid: 2 g/dL	Albumin: 2 g/dL
Creatinine: 200 mg/dL	Oxalic Acid: 60 mg/dL
Hemoglobin: 5 g/L	Bilirubin: 20 mg/dL
Triglyceride:3.89 mmol/L	Cholesterol: 6 mmol/L

None of the substances at the concentration tested interfered in the assay.

【BIBLIOGRAPHY】

1. 1.Jonsson CB，Figueiredo LT，Vapalahti O. A global perspective on hantavirus ecology，epidemiology，and disease［J］. Clin Microbiol Rev，2010，23（2）:412-441.
2. Dzagurova T，Tkachenko E，Slonova R，et al. Antigenic relationships of hantavirus strains analysed by monoclonal antibodies［J］. Arch Virol，1995，140（10）：1763-1773.
3. Kucinskaite-Kodze I，Petraityte-Burneikiene R，Zvirbliene A，et al. Characterization of monoclonal antibodies against hantavirus nucleocapsid protein and their use for immunohistochemistry on rodent and human samples［J］. Arch Virol， 2011， 156（3）:443-456.
4. Hujakka H，Koistinen V，Kuronen I，et al. Diagnostic rapid tests for acute hantavirus infections：specific tests for Hantaan，Dobrava and Puumala viruses versus a hantavirus combination test［J］. J Virol Methods，2003，108（1）：117-122.
5. Araki K，Yoshimatsu K，Ogino M，et al. Truncated hantavirus nucleocapsid proteins for serotyping Hantaan， Seoul， and Dobrava hantavirus infections［J］. J Clin Microbiol，2001， 39（7）:2397-2404.
6. Nichol ST，Spiropoulou CF，Morzunov S，et al. Genetic identification of a hantavirus associated with an outbreak of acute respiratory illness［J］.Science，1993，262（5135）:914-917.
7. Hussein IT，Haseeb A，Haque A，et al. Recent advances in hantavirus molecular biology and disease［J］. Adv Appl Microbiol，2011，74：35-75.
8. Padula, P. J., Edelstein, A., Miguel, S. D., Lopez, N. M., Rossi, C. M., and Rabinovich, R. D. (1998). Hantavirus pulmonary syndrome outbreak in Argentina: Molecular evidence for person-to-person transmission of Andes virus. Virology 241, 323–330.

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IVD	*In vitro* diagnostic medical device
2°C 30°C	Temperature limit
	Do not use if package is damaged and consult instructions for use
EC REP	Authorized representative in the European Community
REF	Catalogue number
Σ	Contains sufficient for <n> tests
	Use-by date
LOT	Batch code
	Manufacturer
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Caution

Hangzhou AllTest Biotech Co.,Ltd.
#550,Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou, 310018 P.R. China
Web: www.alltests.com.cn Email: info@alltests.com.cn

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【BIBLIOGRAFIJA】

1. Jonsson CB, Figueiredo LT, Vapalahti O. A global perspective on hantavirus ecology, epidemiology, and disease [J]. Clin Microbiol Rev, 2010, 23 (2): 412–441.
2. Dzagurova T, Tkachenko E, Slonova R, et al. Antigenic relationships of hantavirus strains analysed by monoclonal antibodies [J]. Arch Virol, 1995, 140 (10): 1763–1773.
3. Kucinskaite-Kodze I, Petraityte-Burneikiene R, Zvirbliene A, et al. Characterization of monoclonal antibodies against hantavirus nucleocapsid protein and their use for immunohistochemistry on rodent and human samples [J]. Arch Virol, 2011, 156 (3): 443–456.
4. Hujakka H, Koistinen V, Kuronen I, et al. Diagnostic rapid tests for acute hantavirus infections: specific tests for Hantaan, Dobrava and Puumala viruses versus a hantavirus combination test [J]. J Virol Methods, 2003, 108 (1): 117–122.
5. Araki K, Yoshimatsu K, Ogino M, et al. Truncated hantavirus nucleocapsid proteins for serotyping Hantaan, Seoul, and Dobrava hantavirus infections [J]. J Clin Microbiol, 2001, 39 (7): 2397–2404.
6. Nichol ST, Spiropoulou CF, Morzunov S, et al. Genetic identification of a hantavirus associated with an outbreak of acute respiratory illness [J]. Science, 1993, 262 (5135): 914–917.
7. Hussein IT, Haseeb A, Haque A, et al. Recent advances in hantavirus molecular biology and disease [J]. Adv Appl Microbiol, 2011, 74: 35–75.
8. Padula, P. J., Edelstein, A., Miguel, S. D., Lopez, N. M., Rossi, C. M., and Rabinovich, R. D. (1998). Hantavirus pulmonary syndrome outbreak in Argentina: Molecular evidence for person-to-person transmission of Andes virus. Virology 241, 323–330.

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
【BIBLIOGRĀFIJA】

1. Jonsson CB, Figueiredo LT, Vapalahti O. A global perspective on hantavirus ecology, epidemiology, and disease [J]. Clin Microbiol Rev, 2010, 23 (2): 412–441.
2. Dzagurova T, Tkachenko E, Slonova R, et al. Antigenic relationships of hantavirus strains analysed by monoclonal antibodies [J]. Arch Virol, 1995, 140 (10): 1763–1773.
3. Kucinskaite-Kodze I, Petraityte-Burneikiene R, Zvirbliene A, et al. Characterization of monoclonal antibodies against hantavirus nucleocapsid protein and their use for immunohistochemistry on rodent and human samples [J]. Arch Virol, 2011, 156 (3): 443–456.
4. Hujakka H, Koistinen V, Kuronen I, et al. Diagnostic rapid tests for acute hantavirus infections: specific tests for Hantaan, Dobrava and Puumala viruses versus a hantavirus combination test [J]. J Virol Methods, 2003, 108 (1): 117–122.
5. Araki K, Yoshimatsu K, Ogino M, et al. Truncated hantavirus nucleocapsid proteins for serotyping Hantaan, Seoul, and Dobrava hantavirus infections [J]. J Clin Microbiol, 2001, 39 (7): 2397–2404.

6. Nichol ST, Spiropoulou CF, Morzunov S, et al. Genetic identification of a hantavirus associated with an outbreak of acute respiratory illness [J]. Science, 1993, 262 (5135): 914–917.
7. Hussein IT, Haseeb A, Haque A, et al. Recent advances in hantavirus molecular biology and disease [J]. Adv Appl Microbiol, 2011, 74: 35–75.
8. Padula, P. J., Edelstein, A., Miguel, S. D., Lopez, N. M., Rossi, C. M., and Rabinovich, R. D. (1998). Hantavirus pulmonary syndrome outbreak in Argentina: Molecular evidence for person-to-person transmission of Andes virus. Virology 241, 323–330.

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
【KIRJANDUS】

1. Jonsson CB, Figueiredo LT, Vapalahti O. A global perspective on hantavirus ecology, epidemiology, and disease [J]. Clin Microbiol Rev, 2010, 23 (2): 412–441.
2. Dzagurova T, Tkachenko E, Slonova R, et al. Antigenic relationships of hantavirus strains analysed by monoclonal antibodies [J]. Arch Virol, 1995, 140 (10): 1763–1773.
3. Kucinskaite-Kodze I, Petraityte-Burneikiene R, Zvirbliene A, et al. Characterization of monoclonal antibodies against hantavirus nucleocapsid protein and their use for immunohistochemistry on rodent and human samples [J]. Arch Virol, 2011, 156 (3): 443–456.
4. Hujakka H, Koistinen V, Kuronen I, et al. Diagnostic rapid tests for acute hantavirus infections: specific tests for Hantaan, Dobrava and Puumala viruses versus a hantavirus combination test [J]. J Virol Methods, 2003, 108 (1): 117–122.
5. Araki K, Yoshimatsu K, Ogino M, et al. Truncated hantavirus nucleocapsid proteins for serotyping Hantaan, Seoul, and Dobrava hantavirus infections [J]. J Clin Microbiol, 2001, 39 (7): 2397–2404.
6. Nichol ST, Spiropoulou CF, Morzunov S, et al. Genetic identification of a hantavirus associated with an outbreak of acute respiratory illness [J]. Science, 1993, 262 (5135): 914–917.
7. Hussein IT, Haseeb A, Haque A, et al. Recent advances in hantavirus molecular biology and disease [J]. Adv Appl Microbiol, 2011, 74: 35–75.

8. Padula, P. J., Edelstein, A., Miguel, S. D., Lopez, N. M., Rossi, C. M., and Rabinovich, R. D. (1998). Hantavirus pulmonary syndrome outbreak in Argentina: Molecular evidence for person-to-person transmission of Andes virus. Virology 241, 323–330.

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Hangzhou AllTest Biotech Co., Ltd.
MedNet EC-REP GmbH