Rotavirus and Adenovirus Rapid Test (Feces)
Package Insert

REF IRAC-627 | English

A rapid, one step test for the qualitative detection of rotavirus and adenovirus in human feces. For professional in vitro diagnostic use only.

【INTENDED USE】

The Rotavirus and Adenovirus Rapid Test (Feces) is a rapid chromatographic immunoassay for the qualitative detection of rotavirus and adenovirus in human fecal specimens to aid in the diagnosis of rotavirus or adenovirus infection.

【SUMMARY】

Acute diarrhea disease in young children is a major cause of morbidity worldwide and is a leading cause of mortality in developing countries. Rotavirus is the most common agent responsible for acute gastroenteritis, mainly in young children.[1] With hospitalized children suffering from acute enteric disease up to 50% of the analyzed specimen were positive for rotavirus.[2] The viruses replicate in the cell nucleus and tend to be host species specific producing a characteristic cytopathic effect (CPE). Because rotavirus is extremely difficult to culture, it is unusual to use isolation of the virus in diagnosing an infection. Instead, a variety of techniques have been developed to detect rotavirus in feces.

Research has shown that enteric adenoviruses, primarily Ad40 and Ad41, are a leading cause of diarrhea in many of these children, second only to the rotaviruses.[3,4,5,6] These viral pathogens have been isolated throughout the world, and can cause diarrhea in children year round. Infections are most frequently seen in children less than two years of age, but have been found in patients of all ages. Further studies indicate that adenoviruses are associated with 4-15% of all hospitalized cases of viral gastroenteritis.[3,4,5,6] Rapid and accurate diagnosis of gastroenteritis due to adenovirus is helpful in establishing the etiology of gastroenteritis and related patient management. Other diagnostic techniques such as electron microscopy (EM) and nucleic acid hybridization are expensive and labor-intensive. With the self-limiting nature of adenovirus infection, such expensive and labor-intensive tests may not be necessary.

The Rotavirus and Adenovirus Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of rotavirus and adenovirus in human fecal specimen, providing results in 10 minutes. The test utilizes antibody specific for rotavirus and adenovirus to selectively detect rotavirus and adenovirus from human fecal specimens.

【PRINCIPLE】

The Rotavirus and Adenovirus Rapid Test is a qualitative, lateral flow immunoassay for the detection of rotavirus and adenovirus in human fecal specimen.

In this test, the membrane is pre-coated with anti-rotavirus antibody on the T1 test line region of the test and anti-adenovirus antibody on the T2 test line region of the test. During testing, the specimen reacts with the particle coated with anti-rotavirus antibody and anti-adenovirus antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-rotavirus antibody and anti-adenovirus antibody on the membrane and generate a colored line. The presence of these colored lines in test line regions indicates a positive result, while their absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【REAGENTS】

The test contains anti-rotavirus antibody and anti-adenovirus antibody coated particles, anti-rotavirus antibody and anti-adenovirus antibody coated on the membrane.

【PRECAUTIONS】

Please read all the information in this package insert before performing the test.

- For professional use only. For *in vitro* diagnostic use only.
- Do not use after the expiration date. Do not reuse the test.
- Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations.
- Humidity and temperature may adversely affect results.
- Wash hands thoroughly before and after handling.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
- Use the test only once. **Keep the test upright while testing. Do not move or turn the test upside down.**
- The kit must not be frozen or used after the expiration date printed on the package.
- Components provided in the kit are approved for use in the Rotavirus and Adenovirus Rapid Test. Do not use any other commercial kit component.

【STORAGE AND STABILITY】

The kit can be stored at room temperature or refrigerated (2-30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

NOTE: It is suggested to use the test within one hour after removing it from the foil pouch.

【SPECIMEN COLLECTION AND PREPARATION】

- The fecal specimen must be collected in clean, dry, waterproof container containing no detergents, preservatives or transport media.
- No dietary restrictions are necessary before using the Rotavirus and Adenovirus Rapid Test.
- Bring the necessary reagents and feces specimen to room temperature before use.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

【MATERIALS PROVIDED】

- Test Cups (with dilution buffer)
- Package Insert
- Droppers

【MATERIALS REQUIRED BUT NOT PROVIDED】

- Timer
- Specimen Container
- Centrifuge
- Pipette

【DIRECTIONS FOR USE】

Before performing the test, allow the test, specimen to reach room temperature (15-30°C), stool samples must be collected following the instruction below.

1. Wash your hands with soap and rinse with clear water.
2. **Collect fecal specimens:**

Collect sufficient quantity of feces (1-2 mL or 1-2 g) in a clean, dry specimen collection container. Best results will be obtained if the assay is performed within 6 hours after collection. Specimen collected may be stored for 3 days at 2-8°C if not tested within 6 hours. For long term storage, specimens should be kept below -20°C.

3. Bring the pouch to room temperature before opening it. Remove the test cup from the foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.
4. **Process fecal specimens:**

- **For Solid Specimens:**

Unscrew the cap of the test cup and take out the specimen collection applicator.

Randomly stab the specimen collection applicator into the fecal specimen in **at least 3 different sites** to collect approximately 50 mg of feces (equivalent to 1/4 of a pea). Do not scoop the fecal specimen.

- **For Liquid Specimens:**


Hold the dropper vertically, aspirate fecal specimens, and then transfer **2 drops of the liquid** specimen (approximately 50 µL) into the test cup containing the dilution buffer.

5. Insert the specimen collection applicator back into the test cup and tighten the cap.
6. Shake the test cup for about **10-15 seconds** to mix well. Leave the cup for reaction for 2 minutes.
7. Remove the plastic limit strip of the test cup.
8. Put the test cup on a clean and level surface, press the cup body from the top to the bottom and start the timer.

NOTE: Keep the test cup upright during the test developing. Do not move or turn the test cup upside down.

9. **Read results at 10 minutes**. Do not read results after 20 minutes.

NOTE: If the specimen does not migrate (presence of particles), open a new test cup, repeat step 4 and centrifuge the diluted specimen contained in the test cup with a clean tube. Pipette 1-1.5 mL of supernatant, dispense into the test cup and insert the specimen collection applicator back into the test cup and tighten the cap. Continue from step 7-8 onwards in the above instructions for use and start the timer.

【READING THE RESULTS】

POSITIVE:

Rotavirus Positive: *A colored line appears in the control line region (C) and another colored line appears in the T1 line region.

Adenovirus Positive: *A colored line appears in the control line region (C) and another colored line appears in the T2 line region.


Rotavirus and Adenovirus Positive: *A colored line appears in the control line region (C) and two other colored lines appear in T1 line region and T2 line region respectively.

***NOTE:** The intensity of the color in the test line region (T1/T2) will vary depending on the concentration of rotavirus or adenovirus antigens present in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE: One colored line appears in the control line region (C). No colored line appears in the test line region (T1/T2).

INVALID: Control line (C) fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

A procedural control is included in the test. A colored line appearing in the control line region(C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

【LIMITATIONS】

1. The Rotavirus and Adenovirus Rapid Test is for *in vitro* diagnostic use only. The test should be used for the detection of rotavirus and adenovirus in human fecal specimens only. Neither the quantitative value nor the rate of increase in rotavirus and adenovirus concentration can be determined by this qualitative test.
2. The Rotavirus and Adenovirus Rapid Test will only indicate the presence of rotavirus and adenovirus in the specimen and should not be used as the sole criteria for the conforming rotavirus and adenovirus to be etiological agent for diarrhea.
3. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of rotavirus or adenovirus infection with low concentration of virus particles.
5. Some nonspecific reactions were documented by the newborns: the presence of both adenovirus and rotavirus lines should be interpreted carefully.

【PERFORMANCE CHARACTERISTICS】

Accuracy

The performance of the Rotavirus and Adenovirus Rapid Test has been evaluated the clinical specimens in comparison with the other rapid test method. The results show that the Rotavirus and Adenovirus Rapid Test has high sensitivity and specificity for rotavirus and adenovirus.

Method		Other Rapid Test		Total Results
Rotavirus Rapid Test	Results	Positive	Negative	
	Positive	49	3	52
	Negative	1	97	98
Total Results		50	100	150

Relative Sensitivity: 98.0% (95%CI:*89.4%-99.9%) *Confidence Intervals
Relative Specificity: 97.0% (95%CI:*91.5%-99.4%)
Accuracy: 97.3% (95%CI:*93.3%-99.3%)

Method		Other Rapid Test		Total Results
Adenovirus Rapid Test	Results	Positive	Negative	
	Positive	48	2	50
	Negative	2	98	100
Total Results		50	100	150

Relative Sensitivity: 96.0% (95%CI:*86.3%-99.5%) *Confidence Intervals
Relative Specificity: 98.0% (95%CI:*93.0%-99.8%)
Accuracy: 97.3% (95%CI:*93.3%-99.3%)

Sensitivity

Rotavirus and Adenovirus Rapid Test can detect the level as low as 50 ng/mL for Rotavirus and 25 ng/mL for Adenovirus.

Dose Hook Effect

There is no dose hook effect with the Rotavirus and Adenovirus Rapid Test.

Precision-Repeatability

Precision-repeatability has been determined by using five specimens: negative, rotavirus low, rotavirus middle, adenovirus low and adenovirus middle positive specimens. The study was performed 5 replicates per day for 5 consecutive days by one operator using 1 lot of Rotavirus and Adenovirus Rapid Test. No difference was observed in intra lot.

Precision-Reproducibility

Precision-reproducibility has been determined by using five specimens: negative, rotavirus low, rotavirus middle, adenovirus low and adenovirus middle positive specimens. The study was performed 3 replicates per day for 5 consecutive days in 3 different sites using 3 separate lots of Rotavirus and Adenovirus Rapid Test (one lot per site), and three operators per site. No difference was observed between days, sites, lots and operators.

Cross-Reactivity

Cross reactivity with following organisms has been studied at 1.0 x 10^7 org/mL. The following organisms were found negative when tested with the Rotavirus and Adenovirus Rapid Test.

Staphylococcus aureus	*Proteus mirabilis*	*Neisseria gonorrhea*
Pseudomonas aeruginosa	*Acinetobacter spp*	*Group B Streptococcus*
Enterococcus faecalis	*Salmonella choleraesius*	*Proteus vulgaris*
Group C Streptococcus	*Gardnerella vaginalis*	*Enterococcus faecium*
Klebsiella pneumoniae	*Acinetobacter calcoaceticus*	*Hemophilus influenzae*
Branhamella catarrhalis	*E. coli*	*Neisseria meningitidis*
Candida albicans	*Chlamydia trachomatis*	

Interfering Substances

The following potentially Interfering Substances were added to Rotavirus and Adenovirus negative and positive specimens.

Ascorbic acid: 20 mg/dL	Oxalic acid: 60 mg/dL	Bilirubin: 100 mg/dL
Uric acid: 60 mg/dL	Aspirin: 60 mg/dL	Urea: 2000 mg/dL
Glucose: 2000 mg/dL	Caffeine: 40 mg/dL	Albumin: 2000 mg/dL

None of the substances at the concentration tested interfered in the assay.

【BIBLIOGRAPHY】

1. WILHELMI I, ROMAN E, SANCHEZ-FAUQUIER A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
2. Cukor, G, Perron, DM, Hudson, R and Blacklow, NR. Detection of Rotavirus in Human Stools by Using Monoclonal Antibody. J. Clin. Microl. 1984, 19: 888-892.
3. Wood, D. J. and A. S. Bailey. Detection of Adenovirus Types 40 and 41 in Stool Specimens by Immune Electron Microscopy. Journal of Medical Virology, 1987; 21: 191-199.
4. Nishio, Osamu, M. Ooseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces. Microbiol. Immunol. 1990, 34(10): 871-877.
5. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens. Journal of Clinical Microbiology, June 1989, 27(6): 1155-1158.
6. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis. Am. J. Clin. Pathol. 1994, 101:742-746.

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